Titan 400_{HP} S400T

(US PATENT 7,510,313)

Xenon Illuminator Operator Manual CE

Sunoptic Technologies®
6018 Bowdendale Avenue
Jacksonville, FL 32216 USA

Customer Service: 904 737 7611
Toll Free 877-677-2832

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

TABLE OF CONTENTS

- 1. INTRODUCTION**
- 2. WARNINGS**
- 3. SPECIFICATIONS**
- 4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS**
 - 4.1 FRONT PANEL**
 - 4.2 REAR PANEL**
- 5. INSTALLATION**
 - 5.1 SETTING UP LIGHT SOURCE**
 - 5.2 CONNECTING LIGHT CABLE**
- 6. OPERATION**
 - 6.1 POWERING UP LIGHT SOURCE**
 - 6.2 SHUTTING DOWN THE LIGHT SOURCE**
 - 6.3 LIGHT BRIGHTNESS CONTROL**
 - 6.4 LAMP LIFE HOURS DISPLAY**
- 7. CLEANING**
- 8. DISINFECTION**
 - 8.1 DISINFECTING CONTROL UNIT**
- 9. MAINTENANCE, SERVICING & REPAIR, WARRANTY**
 - 9.1 LAMP REPLACEMENT**
 - 9.2 WARRANTY**
- 10. END OF PRODUCT LIFE**
- 11. TROUBLESHOOTING**
- 12. CHART OF MEDICAL DEVICE SYMBOLS**

1. INTRODUCTION

Congratulations on the purchase of your new Titan 400 *Xenon Illuminator!*

This user-friendly Xenon Illuminator is a high efficiency light source utilizing state-of-the-art illumination technology. It offers a variety of features such as:

- 5600° K daylight brightness for perfect color definition
- Advanced cooling for superior lamp performance
- Compact and light weight
- Turret which adapts to various types of light guides
- Mechanical iris
- Easy lamp replacement
- Lamp life display indication

This Operator Manual will help you to install the device and optimally integrate it with other components of your system. It will also instruct you how to operate the Xenon Illuminator and how to keep it clean. It will give you maintenance and service guidelines as well as recommendations for best performance results.

INTENDED USE

The intended use of this device is to provide light for fiberoptic cables and instruments – providing light for instrumentation via fiberoptic cables for use in surgical fields.

2. WARNINGS/ CAUTION

Caution

Federal law restricts this device to sale by or on the order of a licensed practitioner.

Caution

To prevent fire or electric shock, do not open or expose the light source unit to rain or moisture. Refer all servicing to qualified personnel only.

Caution

Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Caution

To prevent any potential electro-magnetic interference, do not use any kind of cellular phone near the light source.

Caution

This product should be used only with type BF endoscopic instruments which have been certified according to IEC 60601-1 for medical equipment and IEC 60601-2-18 for endoscopic equipment.

This symbol indicates type BF equipment.

Caution

User must not alter this device in any fashion. Doing so voids all warranties and statements of suitability for any purpose.

Caution

All devices connecting to the Illuminator must be classified as medical equipment. Additional information processing equipment connected to the Illuminator, a Medical System and the operator must determine that all equipment complies with the appropriate end-product standards (such as IEC 60950 or IEC 60065 and the Standard for Medical System, IEC 60601-1-1).

Caution

Always set the intensity control to the minimum level and insert the fiberoptic cable into the unit before turning **on** the power. When light is not required at the surgical site the intensity control should be set to the fully dimmed position. If it becomes necessary to remove the fiber optic cable without turning the unit **off** turn the intensity control to the fully dimmed position.

Caution

The fiber optic cable must be a NON-CONDUCTIVE CABLE. It should not have conductive shielding or any other conductive connection between the patient and equipment. Such connection will impair safety of the equipment. It must be rinsed free of soaking/disinfectant solution and dried before plugging into the light source receptacle. Ensure the optical surface is clean before engaging into the light source.

Caution

Use only approved Sunoptics Surgical[®] cable – part numbers: SG-0080, SG-0083, SG-0087, SG-0088, SG-0089, SG-0090, SL-0080, SL-0083, SL-0087, SL-0088, SL-0089, SL-0090.

Caution

When using high intensity light sources at full output, the recommended distance from the headlight module to the patient is not less than 12 inches (30.5 cm). If using less than 12 inches from the patient, the light intensity must be turned down.

Caution User is responsible for determining if interruption of light output will create an unacceptable risk. If this determination is made arrangements should be taken to reduce the risk.

3. SPECIFICATIONS

Item	Specification
Lamp type	Ceramic type 400 watt Xenon
Power	400 Watt
Color temperature	5600° K
Lamp life	1000 hours (typical)
Lamp replacement	Cartridge replacement from the side of the unit
Light guide adapter	Turret type to fit your choice of four: Storz, ACMI, Wolf, Olympus
Brightness control	Mechanical iris control from the front panel
Input voltage	100-120 V AC, 50/60 Hz 220-240 V AC, 50 Hz
Power consumption	720 watt max
Regulatory Approvals	UL60601, EN 60601-1, EN 60601-1-2, CAN/CSA, C22.2, No 601.1-M90
Equipment Class	Class I, BF-type
Mode of Operation	Continuous operation
Water Resistant	Not Protected Equipment, IPX0
Operating Environment Temperature Relative Humidity Air Pressure	+10° to +40° C (50° to 104° F) 30 to 85% 700 to 1060 hPa
Storage Environment Temperature Relative Humidity Air Pressure	-20° to +60° C (-4° to 140° F) 0 to 95% 700 to 1060 hPa
Dimensions	355 (14)W x 130 (5)H x 250 (10)D mm
Weight	12.4 lbs./5.63 kg

4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS

4.1 FRONT PANEL


Figure 1. Light Source Front Panel

No	Name	Function
1	Turret	Accepts light source end fitting of light guide cable
2	Intensity Control	Mechanically controls the light output
3	Lamp switch	Turns the lamp on and off; illuminates blue when turned on.

4.2 REAR PANEL


Figure 2. Light Source Rear Panel

No	Name	Function
1	AC main power switch and input module	Accepts AC Power cord and switches on AC power
2	Circuit breaker	Over current protection.
3	Hour Meter	Allows tracking of total running time on unit and provides a control on individual lamp hours
4	Fan	Fan 12V DC, airflow for cooling of unit
5	Connector	For potential equalization

5. INSTALLATION

5.1 SETTING UP THE LIGHT SOURCE

Place the light source on a stable surface (cart, counter, stand, etc.).

NOTE

Avoid places where the light source may be splashed with liquid.

Absolutely DO NOT use in any environment with explosive or flammable gases.

DO NOT block the venting grids of the light source or impede exhaust air flow. A minimum of 6"/15cm on each side and 12"/30cm at the back is required, 24"/60cm is preferred. The lightsource should not be placed where its exhaust will influence other devices, nor where the heated air from other devices will influence the lightsource.

Make sure the rear power switch is in OFF position.

Connect AC power cord to the power inlet located on the rear panel of the light source. Make sure power cord is completely seated into the power inlet.

CAUTION

Use only cords provided with the light source.

Plug the AC power cord into a wall outlet using the three (3) prong plug supplied with the unit.

CAUTION

To prevent electric shock, connect power cords of peripheral equipment through medical isolation transformers.

NOTE

When using medical isolation transformer, be sure to check the transformer power ratings.

Make sure that the power cord is connected to the main power with three prong plug (USA use UL2601-1 rated isolation transformers and/or power strips only).

CAUTION

When used in clinical or residential areas near radio or TV units, this equipment may be subjected to radio interference. To avoid adverse electromagnetic effects, DO NOT operate this equipment near RF energy equipment.

5.2 CONNECTING THE LIGHT CABLE

Connect the light cable to instrument and then plug the light cable endfitting into the appropriate port of the turret on the front panel.

NOTE Approved cables for use with the TITAN 400 are: SG-0080, SG-0083, SG-0087, SG-0088, SG-0089, SG-0090, SL-0080, SL-0083, SL-0087, SL-0088, SL-0089 AND SL-0090.

6. OPERATION

6.1 POWERING UP THE LIGHT SOURCE

To operate the light source: Attach the light source to the instrument or headlamp.

Turn on the rear power switch. Turn on the lamp switch. The fan will start. The indicator on lamp switch will light. The lamp will ignite.

6.2 LIGHT BRIGHTNESS DISPLAY AND CONTROL

The light brightness adjustment enables the user to obtain a brighter or darker illumination of the object of observation.

Adjust the light intensity by turning the intensity control knob.

6.3 LAMP LIFE HOURS DISPLAY

The hour meter shows the total elapsed running time of the unit.

7. CLEANING

NOTE *Always disconnect the power cord before cleaning the system.*

The control unit can be cleaned with any cleaning agents, which are used for external cleaning of electric equipment, according to instructions given by the manufacturer of the cleaning solution.

Do not allow excessive moisture or liquids to come in direct contact with the control unit.

Do not use cleaning agents that are not permitted for use with plastics, i.e., ammonia, acetone, salty acids (HCl), etc.

Do not allow cleaning agents or liquids to enter the control unit outlets or vents at side, back or bottom of the unit.

8. DISINFECTION

8.1 DISINFECTING THE UNIT

NOTE	<i>Always disconnect the power cord before cleaning the system.</i>
-------------	--

Use any disinfectant agents which are commonly applied while disinfecting surfaces of electric medical equipment. Such disinfectant agents usually arrive in the form of sprays or damp cloths. Follow the instructions given by the manufacturer of the disinfectant solution.

9. MAINTENANCE, SERVICING AND REPAIR, WARRANTY

Performance of preventive maintenance is not essential. We recommend that once every two months that your biomedical or other hospital technical people remove the lamp module from the light source and use compressed air to blow off any accumulated lint. Regular maintenance can, however, contribute to identifying potential problems before they become serious; thus, enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Defective items of equipment are to be serviced and repaired exclusively by persons authorized by the manufacturer. All repair work shall employ original manufacturer's parts only.

9.1 LAMP REPLACEMENT

WARNING

<p><i>Ceramic lamps are at high internal pressure when cold and at operating temperature; therefore, ceramic lamps may unexpectedly rupture resulting in discharge of hot fragments of quartz and/or glass and metal.</i></p> <p><i>Only handle lamps with protective covers in place. Do not handle lamps without their protective covers unless government approved safety glasses, facemasks (with neck protector, chest protector and gauntlets) are worn.</i></p>
--

CAUTION

<p><i>Always disconnect power cord and turn main switch off before lamp replacement.</i></p>

CAUTION

<p><i>Take utmost care not to hit gas-seal edge of lamp; it may cause gas lamp failure.</i></p>
--

 CAUTION	<i>Do not touch the lamp and patient simultaneously when changing the lamp.</i>
---	--

Use only p/n SSX0040 lamp module for replacement. In order for this light source to perform as designed, the SSX0040 must be used. Failure to use this lamp could void the warranty.

Open lamp compartment door and pull out lamp cartridge.
 Replace lamp cartridge with new unit and log the serial number and hours on the hour meter. Close the lamp compartment door. Re-connect the power cord and turn the light source on according to 6.1. Reset hour meter to zero by depressing hour meter reset button on rear panel.

NOTE	<i>If the tamper seal on the lightsource is broken, the WARRANTY WILL BE VOIDED. There are no other user serviceable parts in the unit.</i>
-------------	--

9.2 LIMITED WARRANTY

Your fiberoptic product(s) carry a 3-year warranty from the date of shipment on workmanship and all defects of material, excluding replaceable sweatbands, headbands, fiberoptic cable and lamps. Should your product prove to have such defects within three years of the shipment, SUNOPTIC TECHNOLOGIES® will repair or replace the product or component part without charge. Should your fiberoptic product(s) need servicing under this warranty, please contact SUNOPTIC TECHNOLOGIES® for return authorization documentation. You should carefully pack unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from SUNOPTIC TECHNOLOGIES®. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

POST WARRANTY REPAIRS: You may return your product(s) for repair, shipping prepaid to the factory. Your product will be inspected and an estimate of repair charges will be submitted to you for approval.

- In the US call: 877 677-2832 (toll free)
- FAX number: 904 733-4832

10. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

In the US a list of recyclers in your area can be found at: <http://www.eiae.org/>.

Please contact customer service to issue a return authorization to return product to manufacturer at the end of product life.

11. TROUBLESHOOTING

<u>Problem</u>	<u>Solution</u>
The lamp switch (refer to 4.1) is not lit.	A. Check that the AC power cord is properly connected. B. Check the circuit breakers. If necessary, reset. C. Power input switch is in the off position
The lamp switch is lit, but Xenon lamp will not ignite.	A. Check to make sure lamp cartridge is seated properly. B. Check that the lamp door is secured. C. Check hour meter; if lamp hours exceed the rated lamp life, replace the lamp (refer to 9.1).

12. CHART OF MEDICAL DEVICE SYMBOLS USED

	Manufacturer
	Date of manufacture (YYYYMM or YYYY)
	“Authorized Representative” in the European Community
	Caution, consult accompanying documents
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner
	Consult Instructions for Use
	CE mark
	Type BF Applied Part
	Not dispose in general waste,
	Caution: Hot Surface
	Caution: Dangerous Voltage
	Product Safety Mark

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